


Lipid status



Sample material: capillary blood / whole blood

Lipid status: Cholesterol, HDL-, LDL-Cholesterol, Triglycerides

1




Blood withdraw with end-to-end capillary:

- HDL 321: 60 μ L

Blood withdraw with ringmark capillary:

- TRI 142: 10 μ L
- CHO 142: 10 μ L

2



Insert 60 μ l end-to-end capillary into reaction tube „R“ and mix well

Wait 5 min.

3



Insert capillaries filled with 10 μ l blood, into TRI 142- and CHO 142-cuvettes, respectively

Wash out blood with Micropipetter

Screw on again turquoise-coloured cap and mix

4



Switch on photometer

Press ON/ENTER

Select TRI and confirm with ON/ENTER

Serial measurement up to 6 samples possible


Lipid status



Sample material: capillary blood / whole blood

Lipid status: Cholesterol, HDL-, LDL-Cholesterol, Triglycerides

5



Measure triglycerides value as per instructions

Photometer saves result

6



Select HDL/CHO and confirm with ON/ENTER

7



Photometer demands for measuring CHO 142 first

Confirm with ON/ENTER

8



Measure cholesterol value as per instructions

Photometer saves result


CHO Cuvette stays into photometer

Lipid status

Sample material: capillary blood / whole blood

Lipid status: Cholesterol, HDL-, LDL-Cholesterol, Triglycerides

9



Continue with HDL-measuring:
Insert reaction tube „R“, with 60 µl capillary (picture 2) into centrifuge

10



Centrifugate 5 minutes

11



Pipette 500 µl supernatant into HDL 321-cuvette

12



Screw on pink cap

Lipid status



Sample material: capillary blood / whole blood

Lipid status: Cholesterol, HDL-, LDL-Cholesterol, Triglycerides

13




Turn HDL 321-cuvette upside down several times

14



Wait 5 min.

15



Now remove CHO 142-cuvette

16



Insert non-treated HDL 321-cuvette into photometer (blank value)

Remove cuvette after signal tone

Lipid status

Sample material: capillary blood / whole blood

Lipid status: Cholesterol, HDL-, LDL-Cholesterol, Triglycerides



17




Insert HDL 321 cuvette into photometer (picture 13)

18



HDL-value will appear on display and be saved by photometer

19



Create measurement report of lipid status with triglycerides, cholesterol, HDL-, LDL-cholesterol and LDL/HDL quotient with DiagWin