Best. Nr. HB 500 HB 2500
Inhalt: 500mL 2500mL

Methode
 Cyanmethämoglobin-Methode^{1,2)}

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Order No. HB 500 HB 2500
Contents: 500mL 2500mL

Method
 Cyan methaemoglobin method^{1,2)}

Sample material
 Capillary blood or EDTA blood
 Use capillary blood immediately. Venous blood can be kept for up to 24 hours at +15°C to +25 °C.

Reagent
 Contents / concentrations of ready-to-use reagent solution:
 Potassium hexacyanoferrate (III) 0.6 mmol/L, Potassium cyanide 0.7 mmol/L, Sodium hydrogencarbonate 18 mmol/L

Safety information
 The reagent is in compliance with CLP-Regulation (EC) classified as hazardous mixture. The reagent solution contains indeed very toxic potassium cyanide, which, however, is considered as not being dangerous due to its extreme low concentration (0.006 %).
 If desired a safety data sheet will be provided.³⁾

Storage and shelf life
 The reagent can be kept in a dark place at a temperature between +15°C and +25°C until the expiry date indicated on the packaging.
 Do not freeze the reagent.

Measurement conditions
Measurement devices: Photometer for rectangular or NP round cuvettes
Meas. wavelengths: 546nm or 560nm
Temperature: Room temperature

Measurement range
 1.0 - 25 g/dL (0.6 - 15.5 mmol/L)

- Tips**
- Store safely away from children.
 - When extracting capillary blood, avoid pressing the finger pulp too hard because otherwise the blood to be extracted is thinned-out by tissue fluid.
 - Cuvettes with a turbid or brownish stained reagent solution cannot be used.
 - If measurement is performed later than after three minutes, the cuvettes must be kept well sealed. Under this condition the colour remains stable for several hours.

Working instructions - Rectangular cuvettes

Pipette into test tube:	
	Analysis
Transformation solution	5.0 mL
Blood	0.02 mL
Wash out the capillary with reagent solution. Mix thoroughly and put into rectangular cuvette. Measure at the earliest after 3 min. against dest. water.	

Working instructions - NP round cuvettes

Pipette into NP round cuvette:	
	Analysis
Transformation solution	5.0 mL
Blood	0.02 mL
Wash out the capillary with reagent solution. Mix thoroughly. Measure at the earliest after 3 min. against dest. water.	

Calculation
 Haemoglobin conc. in blood: $c = \text{Factor} \times \text{Abs.}$

Wavelength	Factor	
	F[g/dL]	F[mmol/L]
546 nm	36.8	22.8
560 nm (Rectangular cuvettes)	42.9	26.6
560 nm (NP round cuvettes)	23.8	14.8

Quality assurance
 For quality assurance we recommend our haemoglobin control **HEM QS**, haemolysate for accuracy and precision control for determination of haemoglobin in normal range.

Reference values

	g/dL	mmol/L
Women	12 - 16	7.45 - 9.93
Men	14 - 18	8.69 - 11.2
Newborn	16 - 25	9.93 - 15.5
Babies	10 - 15	6.21 - 9.31
Toddlers	11 - 14	6.83 - 8.69
Children	12 - 16	7.45 - 9.93

Summary^{2,4)}
 The red blood pigment, haemoglobin (abbreviated Hb) is a protein with iron content which is responsible for the transportation of oxygen in the bloodstream. It consists partly of globulin and the prosthetic haemolysis group. In addition to both major fractions (oxygenated haemoglobin and de-oxygenated haemoglobin), further Hb derivates with a changed haemolysis group (COHb, MetHb) or globulin contents which vary from the standard globulin share (HbA1, HbF) are in the blood.

- Indications / diagnostic significance:
- Recognition of anaemia or hyperglobulia
 - Follow-up and therapeutic controls for anaemia or hyperglobulia
 - Monitoring risk groups for iron deficiency (pregnant women, toddlers, blood donors, haemodialysis patients, sportswomen)

Low HB counts which are not within the reference range are classified as belonging to the symptoms of anaemia and also occur after chronic loss of blood, insufficiently fulfilled additional iron requirements, iron use disorders, intoxication and an array of tumour diseases. Typical symptoms include tiredness and drop in performance. Acute massive haemorrhages normally lead to a considerable drop in the Hb count.

Haemoglobin as a whole is one of the classical POCT (point of care testing) parameters and can be determined in doctors' surgeries using small photometers.

Measurement principle
 Amongst the various determination methods, the cyan methaemoglobin method (also referred to as cyan haemoglobin method) has prevailed and achieved the status of an international reference method^{5,6)}. The concentrations of the reagent's operative elements are defined in a DIN standard⁶⁾.

Through lysis of the erythrocyte membrane, haemoglobin is transferred into the reaction solution, oxidised in methaemoglobin with potassium hexacyanoferrate (III), and converted into cyan methaemoglobin with cyanide. The colour intensity of the cyan methaemoglobin is proportional to the haemoglobin concentration in the sample and is measured photometrically.

Performance parameters

Specificity / interferences
 The physiologically active Hb derivatives (COHb, MetHb etc.) are also registered upon determination. Highly lipaemic samples may disturb the test and falsify HB values which are too high.⁴⁾

Inaccuracy
 The reproducibility was checked using human and control samples.

In series [n = 20]	Average [g/dL]	Standard deviation [g/dL]	VK [%]
Sample 1	5.8	0.07	1.1
Sample 2	11.8	0.12	1.0
Sample 3	16.8	0.17	1.0
From day to day [n = 20]	Average [g/dL]	Standard deviation [g/dL]	VK [%]
Sample 1	5.9	0.09	1.5
Sample 2	11.9	0.13	1.1
Sample 3	16.9	0.20	1.2

Analytic sensitiveness
 Lower detection limit: 1.0 g/dL (0,6 mmol/L)

Comparison of methods
 A comparison of the Diaglobal test HB 500/2500 (y) with a commercially available test (x) resulted in the following correlation according to the Passing/Bablok⁷⁾ process:
 $y = 1.041x - 0.19$
 $r = 0.994$

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